

REQUEST FOR QUOTATION (RFQ)

RFQ Reference: RFQPROC# 2023-0762 - Digital X-ray machine - MHD

Date: 16 June 2023

SECTION 1: REQUEST FOR QUOTATION (RFQ) for the Supply and Delivery of Stationary Digital Xray Machine

International Organisation for Migration (IOM) kindly requests your quotation for the provision of goods, works and/or services as detailed in Annex 1 of this RFQ.

This Request for Quotation comprises the following documents:

Section 1: This request letter

Section 2: RFQ Instructions and Data

Annex 1: Schedule of Requirements

Annex 2: Quotation Submission Form

Annex 3: Technical and Financial Offer

When preparing your quotation, please be guided by the RFQ Instructions and Data. Please note that quotations must be submitted using Annex 2: Quotation Submission Form and Annex 3 Technical and Financial Offer, by the method and by the date and time indicated. It is your responsibility to ensure that your quotation is submitted on or before the deadline. Quotations received after the submission deadline, for whatever reason, will not be considered for evaluation.

Thank you and we look forward to receiving your quotations.

Issued by:

Signature:

Name: Adham Raydan Title: Procurement Assistant Date: 16 June 2023



SECTION 2: RFQ INSTRUCTIONS AND DATA

Deadline for the Submission of Quotation	07 – July - 2023 at 4:00 pm If any doubt exists as to the time zone in which the quotation should be submitted, refer to <u>http://www.timeanddate.com/worldclock/.</u>	
Method of Submission	Quotations must be submitted as follows: □ E-tendering ⊠ Email □ Courier / Hand delivery □ Other Click or tap here to enter text.	
	Bid submission address: iomlebanonbids@iom.int	
	 File Format: PDF or Word 	
	 File names must be maximum 60 characters long and must not contain any letter or special character other than from Latin alphabet/keyboard. 	
	 All files must be free of viruses and not corrupted. 	
	 Max. File Size per transmission: 	
	 Mandatory subject of email: RFQPROC# 2023-0762 - Digital X-ray machine - MHD 	
	 Multiple emails must be clearly identified by indicating in the subject line "email no. X of Y", and the final "email no. Y of Y. 	
	 It is recommended that the entire Quotation be consolidated into as few attachments as possible. 	
	 The proposer should receive an email acknowledging email receipt. 	
Cost of preparation of quotation	IOM shall not be responsible for any costs associated with a Supplier's preparation and submission of a quotation, regardless of the outcome or the manner of conducting the selection process.	
Supplier Code of Conduct	All prospective suppliers must read the UN Supplier Code of Conduct and acknowledge that it provides the minimum standards expected of suppliers to the UN. The Code of Conduct, which includes principles on labour, human rights, environment and ethical conduct may be found at: <u>Supplier Code of Conduct</u> (ungm.org).	
Conflict of Interest	UN encourages every prospective Supplier to avoid and prevent conflicts of interest, by disclosing to UN if you, or any of your affiliates or personnel, were involved in the preparation of the requirements, design, specifications, cost estimates, and other information used in this RFQ.	
General Conditions of	Any Purchase Order or contract that will be issued as a result of this RFQ shall be	
Contract	subject to the IOM General Conditions of Contract for provision of goods/services/transportation/medical services available at https://www.iom.int/do-business-us-procurement.	
Eligibility	Bidders shall have the legal capacity to enter into a binding contract with IOM and to deliver in the country, or through an authorized representative.	



	UN MIGRATION
	 Latest and valid government permits, licenses and company registration Submit <u>any of the following</u> whichever is applicable in the country of business of operation Company's Articles of Incorporation, Partnership or Corporation, whichever is applicable, including amendments thereto, if any License to operate Company Profile (including the names of owners, key officers, technical personnel) showing minimum 1 year of operation supplying similar items as per IOM requirements. Submit duly filled out Vendor Information Sheet and Signed Supplier's Code of Conduct.
Currency of Quotation	Quotations shall be quoted in USD
Duties and taxes	The International Organization for Migration is exempt from all direct taxes, except charges for public utility services, and is exempt from customs restrictions, duties, and charges of a similar nature in respect of articles imported or exported for its official use. All quotations shall be submitted net of any direct taxes and any other taxes and duties, unless otherwise specified below: All prices shall: ⊠ be inclusive of VAT and other applicable indirect taxes Note that IOM shall settle the cost of goods and VAT in USD. □ be exclusive of VAT and other applicable indirect taxes
Language of quotation and documentation including	English
catalogues, instructions, and operating manuals	
Documents to be submitted	 Bidders shall include the following documents in their quotation: ☑ Annex 2: Quotation Submission Form duly completed and signed ☑ Annex 3: Technical and Financial Offer duly completed and signed and in accordance with the Schedule of Requirements in Annex 1 ☑ Other: Formal Quotation based on Supplier's format ☑ Product Brochures or catalogues of offered items ☑ Eligibility Documents
Quotation validity period	Quotations shall remain valid for 90 days from the deadline for the Submission of Quotation.
Price variation	No price variation due to escalation, inflation, fluctuation in exchange rates, or any other market factors shall be accepted at any time during the validity of the quotation after the quotation has been received.
Partial Quotes	 ☑ Not permitted □ Permitted /vendor can quote per line item.
Payment Terms	 ☑ 100% within 30 days after receipt of goods, works and/or services and submission of payment documentation. □ Other
Contact Person for	Focal Person: Adham Raydan
correspondence, notifications and	E-mail address: araydan@iom.int copying beyproc@iom.int Attention: Quotations shall not be submitted to this address but to the address
clarifications and	for quotation submission above.
clarifications	for quotation submission deorer



UN MIGRATION		
Clarifications	Requests for clarification from bidders must be submitted from 16 – June - 2023 until 23- June - 2023. Responses to request for clarification will be	
	communicated via email by 27 June 2023.	
Evaluation method	⊠The contract will be awarded to the lowest price substantially compliant offer in	
	reference to the technical specifications submitted	
	Other Click or tap here to enter text.	
Evaluation criteria – Stage 1	$oxed{\boxtimes}$ Full compliance with all requirements as specified in Annex 1	
	⊠ Full acceptance of the General Conditions of IOM Purchase Order	
	⊠Completeness of Offer	
	⊠Others : Availability of In Country Aftersales Services as part of Warranty	
Evaluation criteria – Stage 2	Quality based on DICOM CXR images to be submitted by the bidders who initially	
	passed the Stage 1 of technical evaluation	
Right not to accept any	IOM is not bound to accept any quotation, nor award a contract or Purchase Order	
quotation		
Right to vary requirement at	Not applicable	
time of award		
Type of Contract to be	Purchase Order	
awarded		
Expected date for contract award.	12 July 2023	
Policies and procedures	This RFQ is conducted in accordance with Policies and Procedures of IOM	
UNGM registration	IOM is encouraging all suppliers to register at the United Nations Global Marketplace	
	(UNGM) website at <u>www.ungm.org</u> . The Bidder may still submit a quotation even if	
	not registered with the UNGM, however, if the Bidder is selected for Contract award	
	of USD 100,000 and above, the Bidder is recommended to register on the UNGM	
	prior to contract signature. For vendors who do not have the technical means to	
	register in UNGM, the UNGM has implemented an assisted vendor registration	
	functionality that allows IOM procurement personnel to add local vendors to the	
	UNGM.	



ANNEX 1: SCHEDULE OF REQUIREMENTS

A. Summary of Requirements

This RFQ refers to the provision of Stationary Digital X-ray Machine as further described in this document titled Annex 1: Schedule of Requirements.

Minimum Technical Specifications for Goods:

Parts	Sub items	IOM standard specification
	1. Statio	nary X ray machine
Brand and Model of the x-ray Machine	Brand/Model	Vendor should specify the brand and model of the X ray machine
X-ray Machine	mA,	650 or more
Generator	KVP,	40 - 150 kV
	Power (KW),	50KW or more
	Automatic exposure control	Present
X-ray tube	Focal spot size of the tube,	0.6/1.2mm (dual)
	Tube Material	Tungsten-Rhenium
	Movability	Multirotating tube (capable to take decubitus view on the table), X-ray tube rotation: ± 180°
	Heat capacity, cooling rate and thermal protection	300KHU or more



		IN MIGRATION
Tube stand	Floor /ceiling amounted tube,	Strong; Horizontal travel range ~ 140 cm; Vertical travel range ~150 cm; Rotation of tube around vertical axis: ± 90°; Floor mounted preferred as it needs less space and no need for strong ceiling (Vendor should specify the type, whether it is floor or ceiling mounted)
Bucky Chest	Bucky (Chest Bucky)	Highest vertical and rotation Bucky movement Distance from chin rest to active area of Bucky allows to fully include the chest x-ray image and avoid cutting of the apices Ergonomic, Strong Bucky stand, which can be well fixed on the floor,
Table	X-ray Table	Four-Way Float Top Table (Minimum Simple, two-way table with locking wheel, especially for small X ray rooms, and if the use is predominantly for chest X rays and requested by the mission).
Grid	Grid Included	Grid ratio, 8:1 or more (the higher is the better)
Exposing switch (Console)	Image exposing Console (If separate from the operator's work station)	Display automatic (Digital) with exposing button
Collimator	Automatic/Semiautoma tic/manual, with light beam centering and timer	(Vendor should specify the type) Automatic better followed by semiautomatic and manual respectively



		IN MIGRATION
Electricity	Electricity requirement	Either Three Phase
source	of the x-ray machine	
requirement		Or
requirement		
		If only single phase electricity is available in the mission, the Single phase machine should be with DQE >80%)
Room	Room size and shape	Vendor should specify the room size
requirement	needed, for stationary	requirement of the X ray machine
	(Room size and x-ray	
	arrangement should	
	allow lateral decubitus	
	view)	
	2. 1	he Detector
_	D	Vendor should specify the brand and model of the
The Detector specification	Brand/model	detector
	One/double detector,	
		detector preferred if double for both table and Bucky instead
		detector
	One/double detector,	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table
	One/double detector, Physical weight	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17'', better than 14X17'' Cesium Iodide (CsI), better than Gadolinium
	One/double detector, Physical weight Physical size of detector	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17'', better than 14X17'' Cesium Iodide (CsI), better than Gadolinium oxysulfide (Gadox)
	One/double detector, Physical weight Physical size of detector Detector Material	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium Iodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the
	One/double detector, Physical weight Physical size of detector Detector Material	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17'', better than 14X17'' Cesium Iodide (CsI), better than Gadolinium oxysulfide (Gadox)
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium Iodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better)
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium lodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum efficiency)	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17'', better than 14X17'' Cesium lodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for single phase machine
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum efficiency) Analog to Digital (AD)	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium lodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum efficiency) Analog to Digital (AD) Conversion	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium lodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for single phase machine 12 bits or more
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum efficiency) Analog to Digital (AD)	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17'', better than 14X17'' Cesium lodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for single phase machine
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum efficiency) Analog to Digital (AD) Conversion	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium Iodide (CsI), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for single phase machine 12 bits or more Wired/wireless image transfer –or wireless through Wi-Fi router, the best specification and latest model
	One/double detector, Physical weight Physical size of detector Detector Material Pixel number Pixel size/ pitch DQE (Detective quantum efficiency) Analog to Digital (AD) Conversion Image transfer method	detector preferred if double for both table and Bucky instead of moving the detector between the Bucky and table Less than 5kg (smaller is better) 17X17", better than 14X17" Cesium lodide (Csl), better than Gadolinium oxysulfide (Gadox) More than or equal to 3000 x 3,000 (the higher the better) Less than or equal to 140 μm (the lower the better) 65% or more for three phase machine and >80% for single phase machine 12 bits or more Wired/wireless image transfer –or wireless through Wi-Fi router, the best specification and latest model of router



		IN MIGRATION
Detector	Charging method	1. direct electric charging preferred
Battery		2. if not, external battery charging with tethering
		option
	If external battery	At least 3 batteries
	charger, number of	
	batteries	
	If external battery	after full charging, can take 200 or more CXRs or
	charger, number of chest	lasts for 8 hours of operation
	x-ray it can take after full	
	charge of the battery	
	If direct electric charging,	If the machine has double detector for both chest
	is detector	and table Bucky separate, integrated fixed detector
		is better.
	integrated(internal)/ or	
	movable (external)	But if single detector, external (movable) detector is
		preferred, so that it can be moved between the
		chest Bucky and the table
	3. Computer	(hardware & software)
Computer	Operating system	windows 10, preferably
(hardware &		
software)	Memory	16GB or more
		· · · · · · ·
	Processor	latest, i7 or above
	Image storage	1TB or higher (can store more than 3000 CXR
		images)
	Image manipulating	Specify brand name of software
	software	Specify Statia name of Software
	DICOM compliance	DICOM 3 compatible & IHE compliant, with ability to
	Dicolarcompliance	network to PACS
	Ports	Ports for USB, DVD/CD burning ability
	Polts	
	Radiographer monitor	minimum 17"
	size	
	HIS/RIS	included and enabled
		included and enabled
	Other supports	Supports worklist modality, storage, printing and
		query and retrieve.
	Biodata fields on the	All biodata fields can be accessible on the display
	display	with blodded fields call be decessible of the display
	•	Other Items
	1	· · · ·
Certification	Registration by regulatory	The machine is approved by authorized regulatory
	body	bodies and manufacturer is legally registered, and
		agent/vender is available in the country,
	Certification document	Vendor can submit any of the following certification
		documents:
		1. Preferably if WHO/SRA certification such as FDA
		or EU certification.
		2. If not, Certification from the country where the X-
		ray machine will be installed +/- CE marking
		3. If 1 and 2 are not available, Certification from the



	U	IN MIGRATION
		country of manufacture of the X ray machine and
		confirmation that the X-ray machine is allowed to be
		installed in the country where the x-ray machine will
		be installed.
Accessories	Lead aprons, skirts and	Vendor to specify the type and quantity of lead
Included	gloves	accessories which can be provided
	(4 Aprons for Technicians,	
	3 Aprons for pregnancy, 3	
	normal Aprons for	
	applicants)	
	Lateral cassette holder	Specify if Lateral cassette holder can be provided
		and included in the package (needed for taking
		lateral decubitus views)
	Child Immobilizer	Specify if include in the package
		specify if include in the package
	Other accessories (such	Include in the optional list the offer for any other
	as x-ray printer, CD	available accessories related to digital radiography
	burners and others)	services
After sale	Availability of technical	Vendor should confirm if technical support is
service by the	support locally	available locally in the country where the X ray
vende		machine is going to be installed
	The capacity of the	The vendor should confirm the technicians available
	technicians	locally are able to fix any machine related issues
	The easy availability of	The vendor should confirm the technicians available
	technical support	locally are accessible in regular and emergency bases
	Response time to critical	
	failure (The equipment	Uninterrupted support with response in few hours
	stops working)	preferred, onsite or remote,
	Response time to non	
	critical failure (the	within 3days
	equipment is still working)	within Sudys
Spare part	Availability of spare parts	
availability	in Stock	All spare parts are available locally
	Time needed to get spare	The vendor should specify how long does it take to
	part, when not available	get spare parts if not available locally
	in stock	(within a week)
Warranty	Duration of warranty	Preferably 2 years
wairancy	Duration of warranty	Preferably 2 years
	The parts/actions	The warranty is comprehensive including the X ray
	included in the warranty	machine tube, detector and other high vacuum
		items. (including changing new detector-if the
		detector fails)
	The service and	24/7 on call for troubleshooting, and regular every
	maintenance agreement	3months (calibrations, software updates)
	included in the warranty	
Availability	Service maintenance	Vender to specify if they have service and
of Service	and insurance coverage	maintenance service, which can be used after
and	services is available	the warranty period, and provide the cost in the
maintenance	after the warranty	optional offer.
services	•	
services after the	period	
	-	



	U	IN MIGRATION
warranty period		
Global reputability of the brand in the sector including the years of experience in the market and field-proven products	Good/Bad/No experience	at least 25 years experience of the manufacturer of X ray machine with no known technical issues
IOM experience on the machine	Good/Bad/No experience	Vendor to specify if similar items have been provided to IOM offices or clinics. Note that IOM internal experience in using the Machine in one of IOM health assessment clinics worldwide will be taken in to consideration.
Image quality checks	DICOM CXR images can be provided for image quality checking	The vender agrees to send at least 5 chest X ray DICOM (Not JPEG)images taken by the same machine brand and model installed at another institution has to be checked by a radiologist to ensure the acceptability of the image it produces.
Users opinion	The vendor has installed the machine to other users, and can provide contact details to collect users' opinion	Vendor to provide contact list of users from where the same brand and model of machine is installed to get users opinion Vendor to provide 5 recommendation letters from centers where the same brand and model of machines are installed (this will be done if IOM Doesn't have direct experience in using the machine if the technical specifications are met.)
Delivery time and place,	Duration of time needed for delivering the machine and certainty, and the delivery is to the site where the machine will be installed,	Please provide delivery lead time



		IN MIGRATION
Installation	Included in the Package	The installation will be done onsite by the
		vendor as part of the package
Users	Included in the package	Training will be given to the users by the vendor
training		onsite, and users guideline will be provided, as
_		part of the package
Assisitng the	Guide in completing the	The vendor can assist to get approval for license
process for	requirements, when	to operate x-ray from Government where the X
getting	needed	ray machine is going to be installed, if needed
license to	needed	Tay machine is going to be installed, if fleeded
operate x-		
ray from		
Gov. if		
needed		
Cost per unit	inclusive of all in the	Amount for each part and the whole package
	package (including	
	tax/VAT)	
	5	. CD Burner
	-	
CD burner	DICOM Compatibility	DICOM part 10 compliant and ability to network
		with PACS
	Automation	Automatically retrieve patient studies from
		networked PACS
	Label Creation	Create customized label with patient biodata
		information using the inbuilt printer. Printing using
		Inkjet or Thermal transfer technology
	Speed	Capable to produce 60 CDs in an hour. Preferably
	2100111	separate bin for input and output of CDs
	DICOM Viewer	Standard DICOM viewer included in the DICOM disc
	Dice Deserviers	created
	Disc Recorders	Single or Dual Disc Recorders. Can choose dual disc
	Software	recorder if the volume is very high Easy and intuitive for the user to use
	SUILWAIE	Easy and intuitive for the user to use
	Media	CD or DVD
	6.)	X-Ray Printer
X-Ray Printer	Does it include the x-ray	Vendor to specify if X ray printer price is included in
	printer	the package. If not, please provide the unit price
		separately.



Delivery Requirements

Delivery Requirements			
Delivery date and time	Goods shall be delivered to the place according to the delivery schedule outlined in each Purchase Order.		
Delivery Terms (INCOTERMS 2020)	n/a – delivery to IOM Lebanon should not include IOM involvement in customs clearance.		
Customs clearance (must be linked to INCOTERM	 Not applicable Shall be done by: Name of organisation Supplier/bidder: delivery to IOM Lebanon should not include IOM involvement in customs clearance. Freight Forwarder 		
Exact Address(es) of Delivery Location(s)	IOM Lebanon Ramlet El Bayda, Nicolas Ibrahim Sursock Street, Beirut, Lebanon		
Distribution of shipping documents (if using freight forwarder)	Not applicable		
Packing Requirements	Standard packaging		
Training on Operations and Maintenance	Manuals and Instruction for Use (IFU) should be included in the standard packaging		
Warranty Period	Option 1: 2 years Standard manufacturer's warranty as applicable to the product Option 2: 3 years Standard manufacturer's warranty as applicable to the product Option 3: 5 years Standard manufacturer's warranty as applicable to the product		
After-sales service and local service support requirements	Must be available in Lebanon Please indicate in your offer the name of the After Sales Service Center if different from Bidder.		
Preferred Mode of Transport	By land		
Other information	This is open to local vendors in Lebanon or vendors with capacity to deliver to IOM without involving IOM on customs clearance.		